




HYSITRON[®]
NANOMECHANICAL TEST INSTRUMENTS

Platform Systems

CONFIDENTIAL


Mechanical Testing Solutions

- Hardness – Modulus
- Dynamic Analysis
- Scratch, Wear, Friction
- Surface Roughness, Imaging
- Environmental Control
- Electrical Properties
- Stress/Strain, Strain Rate





Dual-Head TI 950 TriboIndenter® powered by the *performech*®


Why Test at Nanoscale?



Small Feature Position
Accuracy
Carbide Inclusion - Steel



Thin film characterization
Cohesive fracture – Thin ULK film



Depth and Load Control
Glass Fiber Matrix


Hysitron® Instruments




ATI-8800



TI 950
TribolIndenter




TI Premier



PI 87/85 SEM PicoIndenter





PI 95 TEM PicoIndenter



TS 75 TriboScope

Testing Techniques

Nanoindentation




Hardness, Modulus, Dynamic




HYSITRON®

Tribology


Scratch



Adhesion




Wear




SPM Imaging



*Surface
Roughness*



Positioning



Nanoindentation



Applications




Hardness – Reduced Modulus –
Fracture Toughness

SPM imaging
Metal alloy– NiTi



High Throughput Nanoindentation - XPM


20 x 20 XPM Nanoindentation Grid
400 indents in roughly 100 seconds




HYSITRON®

Optical – SPM – XPM Indentation
Ceramic Matrix Composite – SiC/SiC



Nanowear



Topography




Friction Map



Ramped wear test
Thin film – DLC coating

Nanoscratch



SEM image of the nanoscratch



HYSITRON®





Image Scan Size: 20 μm

SPM image of a nanoscratch


Applications

Scratch Resistance – Film
Adhesion – Coefficient of Friction




Lateral force signal during the scratch


SPM Imaging



40 x 40 μm SPM Image
Glass Beads in Polymer Matrix



3D Image
Glass Beads in Polymer Matrix




1 x 1 μm SPM image
Graphene sheet




HYSITRON®


Environmental Control



xSol – Heating, Cooling, Humidity
Heating to 800 C, Cooling to -100 C, Humidity to 95% RH







High Temperature Ductile Transition
Nanoindentation in Silicon




High Temperature High Humidity
Mechanical Properties
Nanoindentation on Cellulose Acetate

Additional Techniques

- *nanoECR*
 - *Mechanical + Electrical Characterization*
- *nanoDMA III*
 - *Viscoelastic Properties, Creep and Depth Profile*
- *Mechanical + Raman*
 - *Chemical Analysis and Residual Stress Measurement*





HYSITRON®
NANOMECHANICAL TEST INSTRUMENTS

Thank You!

CONFIDENTIAL